

Lotus Consulting presents:

Hydrogen Fuel Analyzer

The Hydrogen Fuel Analyzer from Lotus Consulting provides impressive separations of trace impurities in hydrogen fuel samples. The system features two Varian 3800 Gas Chromatographs configured in a master/slave setup where a single sample is loaded simultaneously into both instruments for a complete assessment of impurities in the hydrogen fuel gas by 5 separate detectors. The system meets the exacting requirements for most of the components in specifications listed in the California Code of Regulations, Title 4, Division 9, Chapter 6 (www.cdfa.ca.gov/measurement/pdfs/HydrogenPropText.pdf).

The fully automated system is designed to be controlled from a single Workstation where a single method controls both gas chromatographs and all five detectors. Samples are loaded through a 16-position automated sampler and passed on to separate valving schemes for each of the group measurements.

For Total Sulfur, the mandated level of 4 ppb V/V requires a large sample volume be cryo-focused and then directed to a column for separation and then to a sensitive and selective pulsed flame photometric detector. The detected sulfurs are then mathematically summed to yield the total result.



Satellite Gas Chromatograph not pictured.

| H₂ Specification | Value |
|--|--------------|
| Hydrogen Fuel Index (minimum, %) | 99.99 |
| Total Trace Gases (maximum, ppm V/V) | 100 |
| Water (maximum, ppm V/V) | 5 |
| Total Hydrocarbons (maximum, ppm V/V) ⁱ | 2 |
| Oxygen (maximum, ppm V/V) | 5 |
| Helium (maximum, ppm V/V) | 100 |
| Nitrogen and Argon (maximum, ppm V/V) | 100 |
| Carbon Dioxide (maximum, ppm V/V) | 2 |
| Carbon Monoxide (maximum, ppm V/V) | 0.2 |
| Total Sulfur (maximum, ppm V/V) | 0.004 |
| Formaldehyde (maximum, ppm V/V) | 0.01 |
| Formic Acid (maximum, ppm V/V) | 0.2 |
| Ammonia (maximum, ppm V/V) | 0.1 |
| Total Halogenated Compounds (maximum, ppm V/V) | 0.05 |
| Particulates Size (maximum, μm) | 10 |
| Particulate Concentration (maximum, μg/L @ NTP) | 1 |

Total Hydrocarbons are readily measured by measuring an aliquot of the sample and directing it to a flame ionization detector without any separation. This detector is a near perfect carbon counter and results are reported as “ppm as Methane”.

Determinations of Oxygen and Argon are problematic in that they do not separate under most chromatographic conditions, making their individual measurement problematic. This analyzer employs a high performance mass spectrometer to determine each analyte based on their different ions (O₂ - 32m/z and Ar - 40 m/z) generated by the mass spectrometer.

Nitrogen, at the levels required, is measurable with either a pulsed discharge detector or by the mass spectrometer (N₂ - 28 m/z) only if it is separated from Carbon Monoxide (CO - 28 m/z). A typical thermal conductivity detector can often measure down to this specification level; however, as these levels are action levels, the system must be capable of measuring much lower to ensure a proper recording of the actual result.

Helium can only be measured with a thermal conductivity detector with nitrogen as the carrier gas. Care is required to ensure that helium elutes prior to hydrogen to avoid inaccuracies yielded when hydrogen elutes first. The mandated level is quite close to the limits for many detectors; only a top performing thermal conductivity detector can achieve accurate results at this level.

Carbon Monoxide and Carbon Dioxide are readily measured at the specification levels with a reduction catalyst after chromatography to convert them to Methane and detection by Flame Ionization Detector.

No detector is presently available to measure Total Halogenates without first separating and then quantitating each halogenate with a mass spectrometer and summing up the group to report as a total. Although this separation process is very tedious, it does provide an accurate assessment of halogenates present.

Water, Formaldehyde, Formic Acid, Ammonia, and Particulate Size are not measurable by gas chromatography at these levels and must be determined by other techniques.

Measurement of these gases requires a very intricate gas chromatograph with multiple valves and columns, five detectors and multiple ovens for columns and valves. Two Varian 3800 Gas Chromatographs paired as a master/slave can accommodate all the hardware needed to yield results from a single sample loading. All valve actuations, temperature control for all valves, columns and detectors, data collection from all five detectors and report generation are set up through a single Workstation and method.

System Operations

Samples and standards are attached to the 16-position automated sampler. Each one is simultaneously loaded into two cryogenic traps and into four fixed volume sample loops. Flow through the sample loops is turned off just before injection to allow the loops to come to consistent pressure. Most analyses can be concurrent, except for the two measurements performed by the mass spectrometer. Here the analyses must be sequential, but all within the time frame of a single method.

A single master method controls all operations of all valve actuations, column oven settings (including temperature programming of two independent column ovens), both cryogenic traps, and all parameters for the five required detectors and their auxiliary flows.

For sulfur, the process starts by allowing the sample to flow through to the mass flow controller prior to trapping to flush out the sample line. The volume loaded onto the trap is determined by the sample flow rate and the time interval for trapping. The trap is isolated during heat-up to ensure a sharp injection into the column.


For speciated halogenates, the process is similar to sulfur, except the column and detector are different.

Measurement of selected fixed gases by mass spectrometry can be performed before or after halogenates by the timing of the gas sample valve actuation. The columns used here are housed in a separate oven for independent temperature setting.

Other fixed gases are measured by injection of a fixed volume sample loop into a column combination selected to strip off heavy interferants and then speciate these analytes with monitoring by two detectors and reduction catalyst all in series.

Total Hydrocarbons are measured by simply directing a fixed volume aliquot directly to the flame ionization detector. As this is executed very rapidly, it can be performed in conjunction with the other measurement involving the flame ionization detector.

Helium can only be determined with a thermal conductivity detector with nitrogen as its carrier and is separated from hydrogen by a long molecular sieve column. Again, potential interferants are stripped away.



System Specifications

Concentrator Traps

- Temperature range: -196 °C to 400 °C
- Maximum heating rate: >300 °C/minute
- Maximum cooling rate: typically >400 °C/minute
- Temperature stability: < 2 °C after 1 minute stabilization
- Temperature overshoot: max. <10 °C, typically <5 °C
- Trap Cryogen usage: < 4 liters per sample
- Trap internal volume: ~100 microliters – cryofocus
~600 microliters – adsorbent trap
- All trap settings controlled/monitored through GC with platinum probe (RTD) and proportional controller (PID)
- Programmable in 5 temperature steps with holds

Automated Sampler

- Standard: 16-position
- Micro-electric actuation, self-aligning
- Independently controlled valve oven
- Maximum temperature limit: 225 °C
- Sample position selected through workstation's sample list
- Position documented in final report and archived with data
- Sample lines heated through control of system

Valving

- Fully automated under time-programmable control of GC
- Total of 14 valves independently actuated
- Valves mounted in heated enclosures
- Micro-electric actuation, easy realignment
- Valco CWE Valves; maximum temperature: 225 °C; some with purged housing and special leak testing
- Valves can be turned on/off 21 separate event times within single method

Sampling

- Sample volumes for Halogenates and Sulfurs set with two independent mass flow controllers with volumes user-selectable through workstation from 5 ml to 1600 ml
- All other volumes are set with fixed volume sample loops
- Loaded sample volume independent of canister pressure

Pneumatics

- Column flows for most column systems employ true Electronic Flow Controller (EFC), not pressure control with computed flows
- Temperature-sensitive flow elements maintained at 45 °C
- Flows automatically adjusted for atmospheric pressure

General

- GC keyboard 11 lines and 35 characters/line for ease of programming and monitoring
- Ethernet communications between GC and Workstation
- Line voltage for GC: 120 V, 20 amperes;
for MS: 120V, 15 amperes

Specifications subject to change without notice

Column Ovens

- Four independently controlled column ovens – two temperature programmable and two isothermal
- Programmable temperature range: -99 °C to 450 °C
- Temperature program rate: 0.1 °C/min to 100 °C/min
- Oven cool-down: 400 °C to 50 °C in 4.5 minutes without cryogen
- Programmable in 7 temperature steps with holds
- Coolant timeout to preserve cryogen when system idle
- Negative temperature programming to save coolant during sample loading

Flame Ionization Detector

- Configured for measurement of Total Hydrocarbons, Methane and CO and CO₂ by post-column reduction to Methane
- Detection to < 0.05 ppm V/V Methane
- Automatic flame-out sensing and reignition
- Electronic flow controllers for supply gases

Reduction Catalyst

- Independently controlled oven to 450 °C
- Detection of CO₂ and CO to < 0.05 ppm V/V
- Electronic flow controller for hydrogen flow

Pulsed Flame Photometric Detector

- Configured for sulfur mode
- Detection < 0.03 ppb V/V H₂S with cryotrapping
- Electronic flow controllers for supply gases
- Selectivity of Sulfur/Carbon > 10⁵
- Linear range > 10²

Thermal Conductivity Detector

- Optimized for detection of helium in hydrogen
- Detection < 10 ppm V/V Helium
- Constant mean temperature setting for filaments
- 20X signal amplification
- Four Tungsten-Rhenium filaments in a Wheatstone Bridge
- Electronic flow control of reference gas

Pulsed Discharge Detector

- Configured for measurement of Nitrogen and Neon and confirmation of Oxygen/Argon, CO₂ and CO
- Detection < 1 ppm V/V Nitrogen
- Detector bypass during elution of Hydrogen

Mass Spectrometer

- Optimized for detection of Halogenates, O₂ and Ar
- Detection of Halogenates < 0.1 ppb V/V
- Detection of Oxygen and Argon < 1 ppm V/V
- Quadruple Ion Trap Design
- Mass range: 10 to 1000 u, time programmable
- Scan rate: 10,000 u/second
- Resolution: better than unit mass (with 10% valley)
- Ion gauge included
- Turbomolecular pumping rate: 280 L/sec
- RJ45 communication between MS and Workstation

Lotus Consulting

310/569-0128 Fax 714/898-7461
email ebramstoncook@msn.com



5781 Campo Walk
Long Beach, California 90803