

Data Integrity with Varian Star Workstation

by:

Randall Bramston-Cook

Lotus Consulting

5781 Campo Walk

Long Beach, Ca 90803

310/569-0128

714/898-7461 - Fax

Copyright 2001 Lotus Flower, Inc.

Computers have become mainstays in the analytical laboratory, especially with chromatographic instruments. They have become data collectors, report generators, system controllers and error monitors. These beneficial tools allow final analytical results to be generated from the detector signal, they permit complete control of the experimental conditions and they can monitor overall system performance. Improved productivity, faster sample turnarounds, better accuracies for results, and easier, more foolproof operations are all gains made over old manual approaches.

However, without some criteria set for monitoring the quality and integrity of the final result, improper use of these computers can yield very inaccurate answers. And without adequate safeguards, opportunities can exist for possible mischief, such as altering sample labeling and archived notes, and even loss of raw data files from accidental overwriting.

The US Food and Drug Administration has been very aggressive in establishing controls to maintain integrity of laboratory data. They have outlined requirements in several sections of the *US Code of Federal Regulations (CFR)*, 4-1-00 Edition (available on-line at www.access.gpo.gov/nara/cfr), including Food and Drug Administration Department of Health and Human Services Title 21, Chapter I, Part 11, Subpart B, and Title 21, Chapter I, Part 58, Subpart G. A summation is provided in CFR, Title 21, Part 11 Subpart B, § 11.10 –

“Persons who use closed systems [environment in which system access is controlled by persons who are responsible for the content of electronic records] to create, modify, maintain, or transmit electronic records shall employ procedures and controls designed to ensure... that the signer cannot readily repudiate the signed record as not genuine.”

Computers can generate results very quickly and very massively; many numbers can be generated that can overwhelm the abilities of the operator to make assure their accuracies and pertinence. Employing a workstation that possesses built-in quality checks can quickly validate this massive data. And ensuring the proper allocations of areas with easy-to-understand algorithms eliminates a tedious manual review of chromatograms and the manual reassignment of areas.

Over the generations of data systems and workstations from Varian (including the CDS111, Vista 401, DS654 and currently Star), Varian has always made strong efforts to protect the integrity of generated results. Varian Star Workstation Version 5.5 is a culmination of those historical efforts and adaptations of new protocols from federal and industry standards, including GLP and ISO9000. The following outlines several of the features in the Varian Star Workstation manifested toward the lofty goal of maintaining data integrity.

This monograph merely lists features available in the Star Workstation and does not include how to accomplish the deeds. The operator is presumed to be knowledgeable enough with the Star Workstation or is willing to explore menus and help screens within the Star Workstation to locate the features from the hints and displays provided.

CFR, Title 21, Part 58 Subpart G, § 58.130 (e) –

“Any change in automated data entries shall be made so as not to obscure the original entry...”

- No Overwriting of Data Files.** Incremented extension XXX is automatically added to yield a unique label if the new file name matches a previous file. The original file cannot be obscured.

Name	Size	Type	Modified
default sample	16KB	Star Workstation Data File	11/9/1999 10:33 AM
default sample001	17KB	Star Workstation Data File	11/9/1999 10:35 AM
default sample002	17KB	Star Workstation Data File	11/9/1999 10:36 AM
default sample003	17KB	Star Workstation Data File	11/9/1999 10:38 AM
default sample004	15KB	Star Workstation Data File	11/9/1999 10:39 AM


- Cannot Change Sample ID.** Once ID is assigned at original data collection, it cannot be altered later and, therefore, cannot be obscured.

```
Title       : Test Run
Run File    : C:\Saturn\MS\GCExamples\Cal_3.run
Method File : C:\star\examples\Is_examp.mth
Sample ID   : 80.0 MP

Injection Date: 12/25/1992 3:58 AM    Calculation Date: 3/26/1998 12:13 PM
```

- Cannot Change Date/Time of Injection.** Once the Injection Date/Time are recorded, they cannot be altered later and, therefore, cannot be obscured.

- Cannot Alter Log Files.** Message, Configuration, Fraction (LC), Run, Revision, Stream Selector (SSV) and Error Logs are automatically generated by Star to document operations within the Workstation. These files can be printed, but cannot be altered later; therefore, they cannot be obscured.






- Cannot Change Sample Notes** - in data file after run is performed; they can be appended with Recalc Notes.


CFR, Title 21, Part 58 Subpart G, § 58.63 (c) –

“Written records shall be maintained of all inspection, maintenance, testing, calibrating, and/or standardizing operations”

6. **All Calibration Data Is Archived With Every Raw Data File.** Details on each standard can be recovered from the data file without an undocumented link to a separate computation method that could have been altered since the original run.




7. **Instrument Maintenance and Repairs Can Be Documented Within Module Notes.** These user-entered notes are saved with every raw data file and cannot be altered later.




8. **Message Log Contains Complete Listing Of System Operations,** including calibration sequence and sample runs. This log can be printed automatically at completion of sequence operations.

	Action	Method	Sample/RecalcList
1	Inject	c:\star\fga.mth	c:\star\fga.smp
2	Inject	c:\star\fga.mth	c:\star\fga.smp
3	Print Message Log		
4			
5			



9. **Injection Counter in Module Notes** provides automatic notification for routine maintenance – for example, septum replacement or column regeneration or column replacement



Injections = 19

10. **Instrument Run Log Documents Instrument Operating Conditions**, even changes made during run – a great assistance in method development. Example run was originally made in 1992, and all GC conditions employed then can be recovered today. Conditions employed to generate data can be verified for compliance with standard operating procedures for the laboratory.


```

3400 GC Event Log:
0.00 Change in Run Status : INJECT
0.00 Initial Detector B Type : FID
0.00 Initial Detector B Attenuation : 8
0.00 Initial Detector B Range : 11
0.00 Initial Detector B Autozero : YES
0.00 Initial Column Temperature (degrees C) : 70
0.00 Initial Injector Temperature (degrees C) : 250
0.00 Initial Detector Temperature (degrees C) : 250
0.00 Initial Auxiliary Temperature (degrees C) : 250
0.01 Initial Pressure A (psi) : 15.5
0.01 Initial Flow A (ml/minute) : 3.7
0.04 AutoSampler Volume (ul) : 2.0
0.04 Start Column Temperature Ramp (degrees C/minute) : 10.0
0.13 Injection Residence Time (minutes) : 0.10
13.01 End Column Temperature Ramp (degrees C) : 200
15.00 Change in Run Status : END
  
```

CFR, Title 21, Part 58 Subpart G, § 58.130 (e) –

“In automated data systems, the individual responsible for direct data input shall be identified at the time of data input”

11. **Automatic Prompt for Operator Name at Start of Data Collection.** Entry for “Operator” is documented on every report in that sequence.



Operator: Randy

```

Method File : C:\My Documents\Brochure
Sample ID   : ELCD Mix
Injection Date: 12/26/2000 1:36 PM   Calculation Date: 1/5/
Operator    : Randy                  Detector Type: ADCB (1
Workstation :                        Bus Address   : 16
Instrument  : PCB Analyzer           Sample Rate  : 10.00 Hz
Channel     : A = ELCD A             Run Time    : 3.198 min
  
```

CFR, Title 21, Part 58 Subpart G, § 58.190 (b) –

“Raw data... may be retained elsewhere provided that the archives have special reference to those other locations.”


12. **StarFinder Documents Location of Archived Files.** When data is transferred off to a storage device, its location can be recorded in StarFinder.






CFR, Title 21, Part 58 Subpart G, § 58.81 (a) –

“A testing facility shall have standard operating procedures in writing setting forth nonclinical laboratory study methods... All deviations...shall be authorized by the study director and shall be documented in the raw data.”

13. **Methods Can Be Set to Require a Password Prior to Saving** – to prohibit unauthorized changes to methods.




14. **To Document Changes in a Method, an Automatic Prompt Can Be Set Up to Add Comments to Revision History** - maintained with the method and with every data file using that method.




CFR, Title 21, Part 11 Subpart B, § 11.10 –

“Procedures and controls shall include...(d) limiting access to authorized individuals.”

15. **Access to Star Security application is password controlled** - limiting access to changes in security procedures and passwords.





16. **Even Specific Applications in Star, such as Interactive Graphics or Standard Reports, Can Be Limited.** When an application is locked, access to it is not permitted. Only after it is unlocked with a password is access allowed.



CFR, Title 21, Part 11 Subpart B, § 11.10 –

“Procedures and controls shall include...(e) use of secure, computer-generated, time-stamped audit trails to independently record time and date of operator entries and actions that create, modify, or delete electronic records.”

17. **Revision Log for Results** - After both checking the box for “Update Revision Log” and enabling the Revision Log in Run Documentation, time and date of operator revisions are independently recorded. This log can be printed as part of every report of results.



```
Revision Log:
12/27/2000 4:59 PM: Calculated results from channel A using method:
'c:\my documents\brochures\data integrity with star workstation\12-26-00 1:37:00 pm elcd mix-b.mth'
12/27/2000 6:20 PM: Calculated results from channel A using method:
'c:\my documents\brochures\data integrity with star workstation\12-26-00 1:37:00 pm elcd mix-b.mth'
12/31/2000 1:40 PM: Calculated results from channel A using method:
'c:\my documents\brochures\data integrity with star workstation\12-26-00 1:37:00 pm elcd mix-b.mth'
```


```
Revision Log:
12/31/2000 2:09 PM: Deleted results calculated on 12/26/2000 9:45 PM
from channel B of ADCB at address 16.
```



CFR, Title 21, Part 11 Subpart B, § 11.10 –

“Procedures and controls shall include... (a) validation of systems to ensure accuracy, reliability, consistent intended performance, and the ability to discern invalid or altered records.”

18. **Validation of System Software Can Be Accomplished Through “Help...About System Control...Validate System Files”**. This test reports that the checksums of the installed applications are identical to those recorded at the factory. Also reported is the software installation history. Separately, the Star Regulatory Compliance Manual, provided on CD, provides detailed procedures for validating computations and report generation.




```
FILE SYSTEM VALIDATION REPORT
File System: Star Workstation 5.51
Workstation: OEMCOMPUTER
Directory: c:\star
Source file: c:\star\star.val
Test Results:
53 files tested
53 files OK
*** File System OK ***
Installation History:
***** Sun Dec 24 10:08:40 2000
Star Chromatography Workstation 5.50
***** Sun Dec 24 10:08:42 2000
Installed: Star Chromatography Workstation 5.50 GC Module Drivers
***** Sun Dec 24 10:08:44 2000
Installed: Advanced Applications 5.50
***** Sun Dec 24 10:10:50 2000
Installed: StarWS 5.50 Update#1 Core Files
Description updated to
Star Workstation 5.51
***** Sun Dec 24 10:10:52 2000
Installed: StarWS 5.50 Update#1 Control Files
***** Sun Dec 24 10:10:52 2000
Installed: StarWS 5.50 Update#1 GC Files
Report Prepared on Sun Dec 31 13:58:46 2000 by _____
Approved: _____
```


Other Features in Star to Assure Data Integrity:

19. **File Names Can Be Assigned Variables** - (mixed with operator-chosen fixed alpha-numeric) for easier file recovery.


Name	Size	Type	Modified
4-5-00 10:34:13 15% co2 20ml fore	30KB	Star Workstation Data File	10/15/2000 9:47...
4-5-00 10:46:48 15% co2 bf@6	37KB	Star Workstation Data File	10/15/2000 9:43...
4-5-00 11:08:14 15% co2 bf@6	34KB	Star Workstation Data File	10/15/2000 9:47...
4-5-00 11:31:06 15% co2 bf@7	38KB	Star Workstation Data File	10/15/2000 9:47...
4-5-00 11:54:00 15%co2 bf@8	33KB	Star Workstation Data File	10/15/2000 9:48...
4-5-00 12:16:51 15% co2 bf@9	37KB	Star Workstation Data File	10/15/2000 9:48...
4-6-00 17:06:26 15% co2 @4.5	31KB	Star Workstation Data File	10/15/2000 9:42...




20. **Files Names Can Be 255 Characters Long**, limited by Microsoft Windows, for clearer labeling and easier recovery of data.

21. **Single Data File Contains Everything Employed To Generate That Data**, including:
 a) raw data for all data channels, b) computed results for all data channels, c) reports for each data channel, d) method used, e) error log, f) run log, g) baseline used in background subtraction, h) snippet of baseline noise monitor, i) calibration data, and j) calibration curves for all data channels. This grouping makes full recovery of data easy and trouble-free and it assures that the original computations can be accurately reenacted later.

22. **Single Method Contains Parameters for Full Instrument Control, Data Collection, Computations, Report Generation, Custom Report, and Printing**. Mistakes in operations are greatly reduced when all parameters are maintained together.



23. **Full Instrument Control from Workstation** – minimizes operating errors when all conditions, including temperatures, flows, AutoSamplers and valves actuations, are combined together in one method.



24. **Auto Start at Specific Time** – For example, the instrument can be set to perform a complete calibration sequence before work shift begins. Calibrations can be performed more often and still keep productivity up.

	Sample Name	Sample Type	Cal. level	Inj.	AutoLink	Rack
1		Autolink			WAIT 3:	
2	Aroclor 1221		1	1	none	1
3	Aroclor 101		1	1	none	1
4	Aroclor 123		1	1	none	1
5	Aroclor 124		1	1	none	1
6	Aroclor 124		1	1	none	1

AutoLink Parameters

Command

WAIT 3:

25. **+1000 Volts to -1000 Volts Electrometer Range on Varian 3800** - No cutoff when signal drifts negatively, or is very large; applicable to FID, TSD, PID, PDD. Range can be set to most sensitive scale and still very large peaks can be measured – over the complete range of most detectors. Eliminates errors when wrong range is chosen.

FID/TSD Detector Full Scale

Front: 1000 V

Middle: 1000 V

Rear: 1000 V

26. **Noise Monitor – Diagnostic for Potential Detector Problems.** Monitor occurs just prior to start of every run. Value is employed in peak detection.

```

Multiplier: 1      Divisor: 1      Unidentified
Baseline Offset: 6248 microVolts
Noise (used): 35 microVolts - monitored before this run
  
```

27. **Peak Detection Based on Meaningful Parameters – Peakwidth and Signal-to-Noise**, not arbitrary setting of “slope sensitivity”. Peakwidths can be obtained from the sample report and S/N is determined from “noise monitor”. Peak areas are computed more accurately and with much less effort when expected peakwidths closely match the actual widths.

Peak Detection

Subtract Blank Baseline

Initial S/N Ratio: 5

Initial Peakwidth: 4 sec

Initial Tangent Height %: 10

Initial S/N Ratio: 5

Initial Peakwidth: 4 sec

Area (counts)	Sep. Code	Width 1/2 (sec)
5273082	BV	5.4
2215177	VV	4.7
2089208	VB	3.2
9577467		

28. **Display of Detector Bunch Points in Graphics** – provides visual indication for adequate number of points to define the peak. Twenty points across the peak normally define the peak properly and yield the best allocation of its area.

Detector Bunch Rate: 8 points (5.0 Hz)

Noise Monitor Length: 64 bunched points (12.8 sec)

Detector Bunch Rate: 8

FID/TSD Detector Full Scale

Front: 1000 V

Middle: 1000 V


Rear: 1000 V

Interactive Graphics




Layout Trace Settings Others Colors

Point


OK Cancel Use Defaults Help




29. **Automatic Tan% Judgement** – Overlapping peaks are problematic when attempting to accurately apportion areas. If the peaks are similar in size, a perpendicular drop to baseline between the peaks is a reasonable estimation. However, when one peak is much smaller than the other, a better guess is to tangent skim. The Tan% parameter allows an automatic judgement based on the relative height of the rider peak to the height of the major peak; less manual intervention in area allocations is then required.



30. **Multiple Reference Peaks** – allow automatic adjustment of peak table for subtle shifts in retention times normally expected in chromatography. If a sizeable change occurs (a shift outside the defined window), a missing reference peak generates an error that can effect a halt in system operations.




31. **Replicate Standard Tolerance** – permits a quality check when multiple standards are employed to generate a calibration curve. If replicates are not consistent (within the defined range), an error is generated that can effect a halt in system operations.




32. **Verification Sample** – After a calibration sequence is performed, one (or more) of the standards can be reexamined as a “verification” sample to provide assurance that the calibration is proper. If a result is outside the tolerance, an error is generated that can effect a halt in system operations.

	Sample Name	Sample Type	Cal. level
1	Control Sample 3	Verification	4
2	Sx 52417281	Analysis	
3			



33. **Calibration Overage** – Sample concentration for an unknown exceeding the highest standard by the tolerance specified will generate an error. This error can effect a halt in system operations.



34. **Error Monitoring** – “Fatal” errors will cause an immediate halt to system operations and a message will be added to the Message Log. “Fatal” errors include: Maximum Non-Fatal Errors Achieved, GC Flameout, Coolant Timeout, Temp Probe Failure, Data Loss Due to Workstation Disconnection, AutoSampler Fault and Gas Saver Timeout. “Non-Fatal errors a) can be documented but ignored, b) can increment an error counter, c) can terminate the sample list and proceed with the next line in the sequence, or d) can halt system operations immediately – options selected by the operator. This group of errors includes: Reference Peak Missing, Internal Standard Peak Missing, Factors Out of Tolerance, Replicate Standard Tolerance, Calibration Range Tolerance, and Verification Sample Tolerance. The operator can choose the maximum number of consecutive faulty runs of the same type that can occur before a fatal error is triggered and operations are halted.

