

Accurate Measure of Injection Volumes with Scion 8400 AutoSampler

Randall Bramston-Cook

Lotus Consulting

5781 Campo Walk, Long Beach, Ca 90803

310/569-0128

randy@lotusinstruments.com

April 30, 2018

Copyright 2008, Lotus Flower, Inc.

Computation of analytical results in gas chromatography usually involves direct comparisons of area counts for target compounds in the unknown sample with corresponding areas for standards, run under identical conditions. Since the two areas are ratioed, accurate aliquots of the injection volumes are not critical, as long as they match. If, for example, the target injection is 1.0 microliters, for both standard and unknown, but instead the injected volume is slightly different, say 1.3 microliters, the change in volume is canceled out in the calculation. Instead, critically important in the experiment is the reproducibility of the injection volume. Any discrepancy in volumes will directly affect the quality of the results. Automated processes for injecting samples into the gas chromatograph greatly improve the precision of the injection volume.

However, circumstances can occur where the accuracy of the specified injection volume plays an important role in the experiment. Examples include:

- Direct comparisons of results with and without internal standard addition.
- Target concentration is bracketed with different injection volumes of the same standard.
- Measurement of extraction efficiencies of purge-trap or solvent extraction with conventional injection of matching liquid standards.
- Multi-point calibration by injecting different volumes of the same standard.

This error occurs because the syringe can have the specified volume in the syringe barrel plus residual sample remaining in the syringe needle just before injection, as illustrated in Figure 1. Typical needles attached to 10 microliter syringes contain 0.7 microliters of sample. When the syringe contents are injected, the actual volume dispensed into the hot injector becomes the desired dosage plus additional amount due to boiling off from residual sample contained in the needle. This amount is very dependent on how long the syringe remains in the injector, the boiling point of the solvent used, and the temperature of the injector. Typically, this residue lost amounts to 0.4 microliters.

To rectify this problem, the Scion 8400 AutoSampler can be preprogrammed to initially insert an air gap into the syringe, a dose of solvent, another air gap, then the sample loading, and finally a third air gap. Figure 3 shows the steps involved. This sequence provides assurance that the full aliquot of the sample is injected, and air gaps and solvent dose flush out the syringe.

After the autosampler sequence is set up to perform pre-injection solvent clean flushes, the syringe needle ends up full of residue solvent, as depicted in Figure 3. This action does not pose any change in the sample dose, as with steps discussed in Figure 1.




Figure 1. Typical sequence of sample loading and injection is shown in this pictorial series. In 1A, the syringe barrel and needle are fully flushed and then the aliquot is set to the desired volume of 1.0 microliter. If the plunger is now pulled back (in 1B), the total volume is shown as 1.7 microliters. After injection, the residue sample is 0.4 microliters, with extra sample pulled out from the heat of the injector. The effective volume injected is 1.3 microliters, or a +30% error from the requested volume.




Figure 2. AutoSampler parameters for operations in Figure 1.




Figure 3. Setting up User Defined parameters for the Scion 8400 AutoSampler gives a sample loading sequence depicted, with insertions of air gaps and solvent flush to provide a full loading of the specified volume. In 2A, the first step is to dose in 1.0 microliters of room air, which becomes a full 1.0 microliter including the volume in the needle. Step 2B then loads up 1.0 microliter of solvent. Step 2C pulls in 1.0 microliter of air to provide a barrier between the solvent and the sample. Step 2D extracts the specified sample volume from the vial. Finally, a third air gap (2E) ensures that no sample persists in the syringe after injection.

GC portal

File Tools Help

GC [10.2.128.5] BrukerGC [10.2.128.3]

SAMPLE DELIVERY MSWS 45X-GC Method Run: 0.00 End: 4.00

CP-8400 Valve Table

Preparation Vial	---	Stroke Count	1
Current Vial	---	Prep. Ahead Time Left	0.00 min
		Syringe Size	10 μ L

Injection Mode: User Defined

Prep. Ahead	<input type="checkbox"/>	Sample Penetration Depth	90 %
Prep. Ahead Delay	0.20 min	Solvent Penetration Depth	90 %

Air Plug After Sample	1.0 μ L	Number of Fill Strokes	0
Sample Air Gap	<input checked="" type="checkbox"/>	Volume for Fill Strokes	1.0 μ L

Viscosity Delay	0.0 s	Fill Plunger Speed	5.0 μ L/s
Pre-Injection Delay	0.0 s	Injection Plunger Speed	5.0 μ L/s
Post-Injection Delay	0.0 s		

Page 1 Page 2 Page 3

Clean Mode

Number of Pre-Injection Solvent Clean Flushes	3
Number of Post-Injection Solvent Clean Flushes	1
Number of Pre-Injection Sample Clean Flushes	0
Clean Solvent Source Vial	I

Internal Standard

Internal Standard	<input type="checkbox"/>	Internal Std. Vial	II
Internal Std. Volume	1.0 μ L	Internal Std. Pause Time	0.0 s
Internal Std. Drawup Speed	5.0 μ L/s	Internal Std. Air Gap	<input checked="" type="checkbox"/>

Page 1 Page 2 Page 3

Solvent Plug

Solvent Plug	<input checked="" type="checkbox"/>	Solvent Drawup Speed	5.0 μ L/s
Vial for Solvent Plug	I	Solvent Pause Time	0.0 s
Solvent Plug Size	1.0 μ L	Solvent Air Gap	<input checked="" type="checkbox"/>

Abort Clean

Clean Vial	I	Number Of Clean Strokes	1
Clean Volume	1.0 μ L	Drawup Speed	5.0 μ L/s

Page 1 Page 2 Page 3

Figure 4. AutoSampler parameters for operations in Figure 3.

The errors realized when using the standard Split/Splitless injection mode is demonstrated in Figure 5. The peak size for lindane is 30 percent larger than when a precise 1.0 microliter is loaded with User Defined parameters listed in Figure 4. The error increases when smaller volumes are specified, as the “hidden” volume in the needle becomes a larger contribution to the effective injection volume.




Figure 5. Comparison of 33 ug/ul lindane chromatograms with Std Split/Splitless Injection Mode versus User Defined Injection Mode applying parameters listed in Figure 4. Targeted injection volume is 1.0 microliter.

A useful capability of properly loading up different injection volumes is the generation of a calibration curve from a single standard. Simply by listing different injection volumes in the sample list from a single vial, achievable results are shown in Figure 6 for lindane. Injection volumes were varied from 0.1 to 0.5 microliters with a 10-microliter syringe.

A similar operation can be performed to bracket sample target concentration with two standards just above and just below the expected concentration simply by stipulating different injection volumes.

Calibration Curve Report

External Standard Analysis
Curve Type: Linear
Origin: Force
y = +5.911398e+006x

Lindane

Resp. Fact. RSD: 8.075%
Coeff. Det.(r²): 0.996915




Figure 6. Calibration curve for 33 µg/µL lindane in iso-octane, generated by using various injection volumes of a single standard with AutoSampler parameters listed in Figure 4.

Lotus Consulting

310/569-0128

Email: randy@lotusinstruments.com



5781 Campo Walk
Long Beach, California 90803