Natural Gas Analysis with Scion 436 Gas Chromatogram

GC ରୁ

One of the most common analysis in gas chromatography has been the measurement of saturated hydrocarbons in natural gas. Analytical results are turned into properties of the gas, including British Thermal Units (BTU), specific gravity, gallons of gasoline per 1000 cubic feet (GPM), and compressibility (Z). A usual protocol often employed for this analysis is GPA Standard 2261-00¹, with fixed gas loop injection, separation by packed columns, and detection with thermal conductivity detector. Calculations follow procedures listed in GPA Standard 2172-09.

The Scion 436 Gas Chromatograph is the perfect instrument for this measurement. It is small and compact, possesses high performance, and can be fully controlled from a remote location. Hardware includes required valving, required columns (configured per method mandates), and thermal conductivity detector. Workstation provides single method for full control of all operating parameters of gas chromatograph, data collection, and report generation of analytical report. Gas chromatograph can be remotely located away from workstation with TCP/IP Ethernet protocol communications.



¹ Gas Processor Association, 6525 East 60th Street, Tulsa, Oklahoma 74145, http://ihsmarkit.com/products/gpa-standards.html

Representative Report:

3 uced Gas ANALY 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	DATE SAMPLED: DATE RECEIVED: GAS TEMP. (°F): GAS PRESSURE (psig): SAMPLED BY: YTICAL RESULTS GROSS HEATING Calculated (Real) Calculated (Real) OT GPM at 14.65	8/22/2018 8/24/2018 60 42 HT G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998 HER	CF) 14.73 DRY 1180.8
<u>uced Gas</u> ANALY 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	DATE RECEIVED: GAS TEMP. (°F): GAS PRESSURE (psig): SAMPLED BY: YTICAL RESULTS GROSS HEATING Calculated (Real) Calculated (Real) OT GPM at 14.65	8/24/2018 60 42 HT G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998	CF) 14.73 DRY 1180.8
MOL % 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	GAS TEMP. (°F): GAS PRESSURE (psig): SAMPLED BY: YTICAL RESULTS GROSS HEATING Calculated (Real) Calculated (Real) Calculated (Real) OT GPM at 14.65	60 42 НТ G VALUE (BTU/SO 14.65 DRY 1160.7 С GRAVITY 0.6998	CF) 14.73 DRY 1180.8
<u>uced Gas</u> ANALY 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	GAS PRESSURE (psig): SAMPLED BY: YTICAL RESULTS GROSS HEATING Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	42 HT G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998	CF) 14.73 DRY 1180.8
Uced Gas ANALY 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	SAMPLED BY: YTICAL RESULTS GROSS HEATING Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	нт G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998 :HER	CF) 14.73 DRY 1180.8
ANALY MOL % 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.20	YTICAL RESULTS GROSS HEATING Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998	CF) 14.73 DRY 1180.8
MOL % 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20 0.30	GROSS HEATING Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	G VALUE (BTU/SO 14.65 DRY 1160.7 C GRAVITY 0.6998	CF) 14.73 DRY 1180.8
MOL % 3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20	Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	DRY 1160.7 C GRAVITY 0.6998 THER	DRY 1180.8
3.11 80.77 1.97 7.25 4.27 0.81 1.05 0.20	Calculated (Real) SPECIFIC Calculated (Real) OT GPM at 14.65	1160.7 <u>C GRAVITY</u> 0.6998 THER	1180.8
80.77 1.97 7.25 4.27 0.81 1.05 0.20	SPECIFIC Calculated (Real) OT	C GRAVITY 0.6998 THER	
1.97 7.25 4.27 0.81 1.05 0.20	SPECIFIC Calculated (Real) OT	C GRAVITY 0.6998 HER	
7.25 4.27 0.81 1.05 0.20	SPECIFIC Calculated (Real) OT	C GRAVITY 0.6998 THER	
4.27 0.81 1.05 0.20	Calculated (Real) OT GPM at 14.65	0.6998 THER	
0.81 1.05 0.20	OT GPM at 14.65	HER	
1.05 0.20	OT GPM at 14.65	HER	
0.20	GPM at 14.65		
0.20	01 11 ut 17.05	18.34	
0.30	Z (dry gas)	0.9968	
0.24			
0.30			
100.00			
A 2166 / ection 1 261-00			
Blauibe isor, tory s			
	0.24 0.30 0.00.00 0.2166 / ection 1 261-00 0.261-00000000000000000000000000000000000	0.24 0.30 0.00.00 2166 / ection 1 261-00 261-00 Sor, ory s the private and exclusive use of Laboratory and its delivery to tion that no representations of warranties expressed or implie information set forth in such report.	0.24 0.30 100.00 2166 / ection 1 261-00 <u>18 could</u> <u>C</u> sor, ory s



Typical chromatogram for Natural Gas, following GPA Standard 2261-00.

Chromatographic Conditions:

Column: Silicone DC 200/500, 30 wt %, 10 m x 1/8 in. OD, plus 0.5 m x 1/8 in. OD plumbed as foreflush and backflush to detector Column temperature: 80.0 °C Flow: 20 ml/min Helium Sample Size: 0.25 ml Detector: TCD Range 0.5 Block Temperature 150.0 °C Filament Temperature: 220.0 °C

Specifications

Scion 436 Gas Chromatograph

Size: Height 57 cm, Width: 32 cm, Depth: 61 cm, Weight: 26.8 kg Environmental Conditions: Operating temperatures: 10 °C to 40 °C Operating humidity (relative): 5% to 95% Line voltage requirements: 230 V ± 10% nominal Column Oven: Dimensions: 23 cm (w) x 11 cm (d) x 28 cm (h) Temperature range: ambient + 4 °C to 450 °C Local Display: TFT full color screen WVGA resolution (800 x 480) Size 23 cm Certifications CSA: C22.2 61010-1 UL 61010-1 IEC: 61010-1 EMC: 47 CFR part 15 **ANSI C63.4** EN 61326 **TCD** Thermal Conductivity Detector Maximum temperature: 450 °C Detectivity: 300 pg/mL (Butane) Linear dynamic range: 10⁶ Operating Mode: constant mean temperature Display of %Balance for troubleshooting Range Selection: 0.05, 0.5 and 5 Valving: Valco valves, actuated with universal electric actuators timeprogrammed through master method, includes operator selection of sample or standard without hardware changes, and sample loop vent to stabilize injection loop pressure. Data Collection, System Control and Communication Ethernet: Protocol: TCP/IP Data rate: 100 Mbps Control: GC control and method parameters Data Processor: Windows 10/32 bit, Intel Core i5 7400 3 GHz, 4 gB memory, 1 tB solid state drive, DVDRW, keyboard, mouse, 27" LCD monitor. Natural Gas Report: Mol%, BTU, Specific Gravity, GPM, Compressibility Factor

Lotus Consulting

310/569-0128 semail randy@lotusinstruments.com



5781 Campo Walk Long Beach, California 90803